

Sports Medicine Grafts

Precision Grafts for Peak Performance

Advanced Human Tendon and
Dermis Solutions in Sports Medicine



Sports Medicine Grafts

Expertly crafted for orthopaedic and sports medicine applications, Tissue Regenix's grafts are processed aseptically using the proprietary BioRinse™ technology, frozen, and sterilised through low-dose, low-temperature irradiation.

Key Features and Benefits

Key features include:



Five-year shelf life



Stored frozen (-40 - -80°C)



High sterility assurance level (SAL) of 10⁻⁶

Product Options

Product Code	Description
311001-TL	Bone-Tendon-Bone (Whole)
311002-TL	Bone-Tendon-Bone (Hemi)
311004-TL	Achilles Tendon w/ Bone Block
311006-TL	Anterior Tibialis Tendon
311007-TL	Posterior Tibialis Tendon
311008-TL	Semitendiosus Tendon
311009-TL	Gracillis Tendon
311010-TL	Peroneus Tendon



Enhance surgical outcomes in sports medicine with our grafts. Reach out for details.

Contact



- info@biostor.eu
- www.biostor.eu





Entering the US regenerative medicine market in 2016, Tissue Regenix now partners with Biostór Ireland, extending their reach into the EU. This partnership ensures Tissue Regenix's innovative regenerative medical technologies are now readily accessible through Import Tissue Establishment, Biostór Ireland, to advance tissue repair and regenerative treatments throughout the European Union.